PULPDENT®

ACTIVA™ PRESTO
BioACTIVE HIGH VISCOSITY • LIGHT CURE
FLOWABLE COMPOSITE

CRYSTA™
MCP TECHNOLOGY
U.S. PATENT NO 11,116,701 B2

XP-VPF-IN-05w
REV: 10/2025

PRODUCT DESCRIPTION

ACTIVA™ PRESTO – BIOACTIVE is an esthetic, light cure material that contains calcium, phosphate, and fluoride in a durable, wear-resistant resin matrix indicated for all restorative procedures. The material is stackable, does not slump, and holds its shape prior to light curing.

ACTIVA PRESTO is tougher than traditional composites. It contains a patented, shock-absorbing, rubberized-resin component that resists chipping and fracture. This increases flexural strength while maintaining the esthetics, hardness, compressive strength, and durability required of dental composites.

ACTIVA PRESTO is a urethane-based resin that contains no Bisphenol A, no Bis-GMA, and no BPA derivatives.

INDICATIONS

Recommended as a filling material for pits, root surface cavities, and Class I, II, III, IV and V restorations where there is no pulpal involvement.

CONTRAINDICATIONS

Not indicated for direct placement on the exposed pulp.

For any potential allergies, refer to associated MSDS.

INSTRUCTIONS FOR USE

1. Isolate and prepare tooth to receive a restoration. Bevel and round off enamel margins. In Class V lesions, bevel or undercut enamel.
2. Place appropriate pulp protection, if indicated.
3. Etch and apply a bonding agent following the manufacturer's instructions for use.
4. **ACTIVA PRESTO** is a dentin and enamel replacement material. Place a 19-gauge applicator tip on the syringe and apply in 2mm increments, light curing for 20 seconds between each layer.
5. Finish and polish as usual.

STORAGE AND HANDLING

• Store tightly sealed in original container at cool room temperature. Avoid direct light, extremes of temperature, contamination and sources of ignition.

• Re-cap immediately after use.

CAUTION

Uncured material may cause eye or skin irritation on contact. Dental professionals should wear safety glasses and surgical gloves. Do not exceed manufacturer's recommended curing time for the light you are using.

REPORTING OF SERIOUS INCIDENTS

Any serious incidents involving the device should be reported to the manufacturer and the Competent Authority of the Member State in which the incident occurred, or other appropriate regulatory body.

Reports of serious incidents should be sent to:

Pulpdent Corporation

80 Oakland Street, Watertown, MA 02472 USA

Tel: 617-926-6666 / Fax: 617-926-6262

email: pulpdent@pulpdent.com

Note: Apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with each patient. For additional information, refer to: <https://www.fda.gov/medical-devices/dental-devices/multi-ple-use-dental-dispenser-devices>

Note: Refer to the Safety Data Sheet for a list of materials present in this product.

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